

Master Blend #5 Non-Acid Bowl Cleaner M2210065

A super concentrated quaternary ammonia detergent and disinfectant. The half gallon containers are designed to be used with the Master Blend Dispensing System and color coded quart bottle system. Master Blend features SafeLink Adaptive Closed-Loop™ Technology. It prevents leaks and is safer for employees and the environment. Economical and easy to use.



WALL MOUNT UNIT
#HYDSP3066



PORTABLE UNIT
#SL1SGAP

Optional Equipment

5 PRODUCT RACK
#ACP627

STANDARD FOAMER
#DEMA52339

PREMIUM FOAMER
#DEMA52336



Where To Use

Athletic facilities	Restrooms	Schools	Nursing Homes
Health clubs	Offices	Restaurants	
Hospitals	Colleges	Institutions	

Dilution

2 oz.. per gallon.

Directions

Master Blend Non Acid Bowl Cleaner #5 is designed to be diluted at 2 oz. per gallon through the Simoniz Master Blend Dilution Systems. The Master Blend systems are preset for a 2 oz. per gallon dilution for filling quarts. A light blue tip must be added to the wall mount unit high flow fill side in order to fill mop buckets at the correct dilution. Simply attach the half gallon to the dispenser via the cap at the end of the draw tube (wall mount) or attach the portable unit directly on to the half gallon and fill either quarts or mop buckets. If using this product with another dispenser, use the appropriate tip for the desired dilution. (Please note, in order to use Master Blend products with other dispensers, part number C38CS4W must be attached on the pickup tube.)

To clean and disinfect toilet bowls - remove or expel over the inner trap residual bowl water prior to applying use solution. Use 1 oz.. of product in toilet bowl. Squeeze around bowl and under rim. Then scrub thoroughly with brush. Let stand 10 minutes then flush.

For urinals - add 1 oz.. of product. Then scrub thoroughly with a brush. Let stand 10 minutes.

For disinfecting - Remove gross food particles and heavy soil deposits, then thoroughly wet surfaces. Apply diluted solution by mop, sponge, cloth or spray on hard nonporous surfaces for a minimum contact time of 10 minutes.

Mildew static instructions - Apply diluted solution by spraying, making sure to wet all surfaces completely. Let air dry. Repeat application weekly or when growth reappears.

HIV-1 (Aids Virus) – The areas to be disinfected must be pre-cleaned of all blood and body fluids. Please read label instructions for use in HIV-1 disinfection.

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

Technical Specifications

COLOR:	Fluorescent pink	EPA Reg. No.:.....	47371-37-18305
FORM:	Liquid	STORAGE(Unopened):	.1 Year @R.T.
SCENT:	Citrus	BIODEGRADABLE:	... Yes
pH:	12.0	WEIGHT/GAL:	8.8 lbs.
FLASH POINT:	None	SHIPPING CLASS:.....	None
VISCOSITY:	Slightly viscous	PACKAGING:	6x64 oz..
FOAM:.....	Low	HOSPITAL GRADE:	Yes
TYPE:	Super Concentrated	HIV CLAIM:	Yes
	Non-Acid Bowl Cleaner		

