

SIMONIZ SIMO-TAB

Disinfectant Tablets
with Enhanced C-diff Kill Time



THE SIMO-TAB C.diff ADVANTAGE

Simo-Tab C.diff deodorizing, disinfecting and sanitizing tablets are the convenient and easy way to kill C.diff and many more organisms **in 4 minutes**. At 6.55 gram per tablet, it takes just 2 tablets per quart to make a powerful disinfectant solution. A single bottle of 120 tablets will yield 60 quarts. *One complete system for infection prevention.*

How It Works:

Simo-Tab C.diff tablets are bleach-free and utilize an active ingredient called sodium troclosene which is much more stable than sodium hypochlorite and releases chlorine at a much slower rate. This means that Simo-Tab C.diff tablets produce a solution that maintains disinfecting power for 7 days in a sealed container and 24 hours in an open container, vs. bleach solutions that become inactive after a day. In solution, sodium troclosene has a pH of around 6.0. This mild pH means that Simo-Tab C.diff tablets will be much less aggressive on surfaces, including floor finishes, than traditional bleach products which regularly have a pH of 11+.

The **Simoniz Electrostatic Sprayer** covers a surface area very quickly and effectively more than traditional manual cleaning that are time-consuming. It coats 360° of all surfaces with 3 times more the coverage than conventional sprayers. It's a simplified, faster and more thorough way to sanitize and disinfect to prevent illness.

Treatable Pathogens:

- Staphylococcus aureus-MRSA & GRSA
- Norovirus
- Influenza H1N1
- C.diff
- Staphylococcus aureus
- Salmonella enterica
- Klebsiella pneumonia
- Staphylococcus epidermidis
- Escherichia coli O157:H7
- Canine Parvovirus
- Canine Distemper Virus
- Feline Calicivirus

Targeted Applications:

- Hospitals
- Clinics
- Restrooms
- Athletic Facilities
- Day Care Centers
- Nursing Homes
- Hotels
- Food Services
- Kennels
- Veterinary Clinics

Simo-Tab C.diff disinfectant tablets are more effective at lower doses than quaternary ammonium compounds, bleach, hydrogen peroxide, and peracetic acids. Germs cannot build up a resistance to Simo-Tab C.diff disinfectant tablets and the solution is not deactivated by organic matter or cloth as other disinfectants can be. Boasting a 3 year shelf life, Simo-Tab C.diff disinfecting and sanitizing tablets ship DOT non-hazardous, are small and convenient and will help you reduce transport space, storage space, and wasteful packaging.



SIMO-TAB

Disinfecting, Deodorizing and Sanitizing Tablets

Packaging:

120 tablets/jar, 2 jars/case

C6000002

SIMONIZ ELECTROSTATIC SPRAYER

Packaging:

1 Unit

ECSPRAYER

Features of Simo-Tab Tablets:

- A SUSTAINABLE and broad spectrum disinfectant, virucide and mildewcide
- EFFECTIVE against C. diff and MRSA in 4 minutes
- DISINFECTS and PROTECTS against odors and odor causing bacteria
- USE with Simoniz Electrostatic Sprayer, spray bottle or mop and bucket
- NEUTRAL pH
- EPA registered 71847-6-18305.

Features of Simoniz Electrostatic Sprayer:

- SPRAYED DROPLETS wrap around all objects, coating 360° of all surfaces 800-1000 sq. ft., **3 times the coverage** of conventional sprayers
- TAKES ADVANTAGE of the negative charge on room surfaces by putting a positive electrical charge on the droplets
- CHARGE allows the droplets to attract to all room surfaces like a magnet
- EASY-TO-USE Battery powered – 4 hour run
- PORTABLE and CORDLESS

