

A concentrated, neutral pH, multipurpose germicidal detergent and deodorant that disinfects, cleans and deodorizes in one labor saving step. Effective against MRSA (Methicillin Resistant Staphylococcus aureus) Trichophyton mentagrophytes, Candida albicans, HIV-1 (AIDS virus), Influenza A/Hong Kong, Herpes Simplex Type 1, Herpes Simplex Type 2, Vaccinia, Rubella (German Measles), Adenovirus Type 4, Canine Distemper and Escherichia Coli (ATCC & antibiotic resistant strain). See efficacy data sheet for additional pathogens.



Where To Use

| | | | |
|-----------|-----------------|-----------------|--------------------|
| Hospitals | Nursing Homes | Supermarkets | Veterinary Clinics |
| Schools | Office Building | Delicatessens | Dumpster Areas |
| Restrooms | Hotels | Finished Floors | Airports |
| Kennels | | | |

Dilution

1/64 (2 oz./gallon)

Directions

General disinfectant - Remove heavy soil deposits from surface, then thoroughly wet surface with a solution of 2 ounces of the concentrate per gallon of water. Apply with cloth, mop sponge, coarse spray or soaking. Let solution remain on surface for a minimum of 10 minutes. Rinse and allow to air dry.

Toilet bowls - Swab bowl with brush to remove heavy soil prior to cleaning and disinfecting. Clean by applying diluted solution around the bowl and up under the rim. Stubborn stains may require brushing. To disinfect, first remove or expel over the inner trap the residual bowl water. Pour in 3 ounces of the dilution. Swab the bowl completely using a scrub or mop, making sure to get under the rim. Let stand for 10 minutes or overnight then flush.

For decontamination against HIV-1 (AIDS virus) of surfaces/objects soiled with blood/body fluid. Disposable latex or vinyl gloves, gowns, face mask or eye coverings, must be worn during all cleaning of body fluids, blood and decontamination procedures. Blood and body fluids must be thoroughly cleaned from surfaces and objects before application of disinfectant. Effective against HIV-1 (AIDS virus) on hard nonporous surfaces in the presence of a moderate amount of organic soil (5% blood serum) (providing 1800 ppm of active quaternary) for a contact time of 4 minutes at room temperature. Use a 10 minute contact time for disinfection against all other bacteria, fungi and viruses claimed (listed on product data sheet).

Technical Specifications

| | |
|--|------------------------------------|
| COLOR: Transparent | EPA Reg No:47371-131-18305 |
| FORM: Liquid | STORAGE(Unopened): 1 Year @R.T. |
| SCENT: Herbal | BIODEGRADABLE: ..Yes |
| pH: 7.5 | WEIGHT/GAL:8.4 lbs. |
| FLASH POINT:.. None | SHIPPING CLASS:None |
| VISCOSITY: Water Thin | PACKAGING:4x1, 5 gal, 55 gal |
| FOAM:..... Low/Medium | HOSPITAL GRADE:....Yes |
| | HIV CLAIM.....Yes |
| TYPE: Herbal Scented Neutral pH Quat | |

